



Salud, Independencia, Alegría



plaza médica[®]
Especialistas en Equipamiento Médico

CATÁLOGO DE SILLA DE RUEDAS

Av. Emancipación 236, Lima / T. (01) 696-8129
WhatsApp: 995 225 764 / 928 396 811 / 930 377 430 / 910 753 085
www.plazamedicaperu.com
ventas@plazamedicaperu.com



plaza médica[®]
Especialistas en Equipamiento Médico



Salud, Independencia, Alegría

Silla de ruedas

Silla de ruedas eléctrica

Cómoda

Silla de ducha

Caminante

Andador

Bastón y muleta

Cama de hospital

Scooter de movilidad y silla de ruedas eléctrica.

silla de ruedas y rampa para scooter



Silla de ruedas y Sco

Silla elevadora



CCW01 
Chromed steel frame
fixed armrest, fixed footrest
solid castor, solid rear wheel

MAIN PARAMETERS (CM)

MODEL NO											CTNS 20' FCL	CTNS 40' HQ FCL
CCW01	60	41	59	20	51	87	100	94x23x89	17.4	19.6	144	372
CCW06	55	35	55	15	47	75	75	88x23x89	16.2	18.4	164	410
CCW07	60	41	59	20	50	87	100	94x25x89	17.3	19.5	132	336
CCW02	60	41	59	20	50	87	100	94x23x89	17.3	19.5	144	372
CCW03	60	41	59	20	50	87	100	94x23x89	17.3	19.5	144	372



CCW06 ▲

Children type, chromed steel frame, fixed armrest
fixed footrest, solid castor, solid rear wheel



CCW07 ▲

Powder coating steel frame, fixed armrest, fixed
footrest, solid castor, solid rear wheel



CCW02 ▲

Chromed steel frame, powder coating cross bar, fixed
armrest, fixed footrest, solid castor, solid rear wheel



CCW03 ▲

Chromed steel frame, fixed armrest, fixed footrest
solid castor, solid rear wheel



CCW94
Powder coating steel frame, double cross bar, fixed inclined armrest, fixed footrest, solid castor, solid rear wheel



CCW09
Chromed steel frame, double cross bar, drop back handle with united brake, fixed armrest, fixed footrest, solid castor, pneumatic MAG rear wheel



CCW05
Chromed steel frame, fixed armrest, fixed footrest, solid castor, solid MAG rear wheel



CCW84
Powder coating steel frame, double cross bar, fixed armrest, fixed footrest, solid castor, solid rear wheel

MAIN PARAMETERS (CM)

MODEL NO.												
	W	D	H	W	H	W	(kg)	(mm)	HW	SW	CTNB 20° FCL	CTNB 10°/90° FCL
CCW84	61	41	59	20	49	87	100	94x23x89	17	19	144	372
CCW05	60	41	61	20	50	88	100	94x25x90	18.7	20.9	140	336
CCW09	60	41	40	15	45	89	100	82x28x72	14.6	16.7	168	387
CCW84	59	46	59	20	47	91	100	93x27x93	16	18	108	272
CCW11	60	41	59	20	47	91	100	94x28x92	17.4	19.6	108	272
CCW103	62	39	20	61	45	84	100	97x31x85	16.1	18.3	88	192
CCW06	61	41	61	20	47	91	100	94x28x92	16.7	20.9	108	272
CCW13	60	41	59	20	47	91	100	94x28x78	17.4	19.6	108	336



CCW11

Chromed steel frame, double cross bar, fixed armrest
fixed footrest, solid castor, solid rear wheel



CCW103

Chromed steel frame, double cross bar, fixed armrest
fixed footrest, pneumatic front wheel, solid rear castor



CCW08

Chromed steel frame, double cross bar, fixed armrest
fixed footrest, solid castor, solid MAG wheel



CCW13

Chromed steel frame, double cross bar, drop back handle
fixed armrest, fixed footrest, solid castor, solid rear wheel



CCW14

Powder coating steel frame, detachable armrest
detachable footrest, solid castor, solid rear wheel



CCW16

Chromed steel frame, detachable armrest, detachable
footrest, PU castor, pneumatic rear wheel



CCW15

Chromed steel frame, detachable desk armrest, elevating
footrest, PU castor, pneumatic rear wheel



CCW17

Powder coating steel frame, detachable armrest
detachable footrest, PU castor, pneumatic MAG wheel

MAIN PARAMETERS (CM)

MODEL NO											
	W	D	H	W	H	H	H	W	H	W	H
CCW14	60	41	59	20	50	90	100	80x28x91	17.5 19.5	126	340
CCW16	60	41	61	20	50	90	100	80x28x91	18 20	128	340
CCW15	60	41	61	20	50	90	100	80x28x91	16.9 19.1	126	340
CCW17	80	41	61	20	50	90	100	80x28x91	20 22	126	340
CCW18	58	41	30	20	50	89	100	71x25x90	14 16.2	180	406
CCW20	62	43	20	20	50	94	100	74x26x78	13 15.2	216	432
CCW19	60	41	59	20	50	87	100	80x24x89	17.6 19.8	172	390
CCW85	60	41	61	20	51	89	100	80x24x91	19.1 21.31	172	360



CCW18

Powder coating steel frame, detachable desk armrest, detachable footrest, solid castor, pneumatic MAG wheel



CCW20

Chromed steel frame, fixed armrest, detachable footrest, drop back handle with united brake, solid wheel



CCW19

Chromed steel frame, fixed armrest, detachable footrest, solid castor, solid rear wheel



CCW85

Powder coating steel frame, fixed armrest, detachable footrest, solid castor, solid MAG wheel



CCW22 ▲

Powder coating steel frame, flip-up armrest, detachable footrest, solid castor, solid MAG wheel



CCW27 ▲

Chromed steel frame, detachable desk armrest, elevating footrest, solid castor, pneumatic MAG wheel



CCW23 ▲

Powder coating steel frame, flip-up desk armrest, detachable footrest, solid castor, solid rear wheel



CCW24 ▲

Chromed steel frame, flip-up desk armrest, detachable footrest, solid castor, quick release pneumatic rear wheel, drop back handle

MAIN PARAMETERS (CM)

MODEL NO.										GTNS-20' FCL	GTNS-42' HD FCL
CCW22	57	41	61	20	50	93	100	82x28x95	19.7 21.9	130	310
CCW27	58	41	40	20	50	91	100	78x28x91	17.2 19.4	180	366
CCW23	56	41	59	20	50	90	100	83x28x95	17.9 20.1	130	310
CCW24	58	41	61	15	50	93	100	80x28x76	16 18	168	402
CCW25	71	56	60	20	54	94	125	83x35x93	22.6 24.8	112	245
CCW26	70	51	60	20	50	92	125	80x28x93	20 22	130	310
CCW104	78	61	30	20	50	96	125	86x37x83	19.3 21.5	144	324
CCW105	80	61	61	20	50	96	125	86x37x83	22.2 24.4	100	240



CCW25

Powder coating steel frame, 22" seat width, flip-up armrest, detachable footrest, rubber foam castor, pneumatic MAG wheel



CCW26

Chromed steel frame, 20" seat width, detachable armrest, detachable footrest, solid castor, pneumatic rear wheel



CCW104

Powder coating steel frame, detachable desk armrest, detachable footrest, 24" seat width, foldable push bar, solid castor, pneumatic MAG wheel



CCW105

Powder coating steel frame, detachable desk armrest, detachable footrest, 24" seat width, foldable push bar, solid castor, pneumatic rear wheel



CCW04

Chromed steel frame, fixed armrest, fixed footrest, solid castor, solid rear wheel, hard seat cushion

MAIN PARAMETERS (CM)

MODEL NO.	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair
CCW04	61	41	59	20	50	89	100	94x31x93	18.8 20.8	108	220	
CCW123	60	41	61	20	50	90	100	80x28x91	19.1 21.2	126	340	
CCW10	60	41	59	20	51	91	110	94x28x93	20 22	108	220	
CCW21	60	41	61	15	51	89	120	94x28x89	19 21	120	248	



CCW123

Powder coating steel frame, detachable armrest
detachable footrest, solid castor, solid MAG
rear wheel



CCW10

Chromed steel frame, fixed armrest,
Fixed footrest, solid castor, solid rear wheel,
Hard cushion, double cross bar



CCW21










Chromed steel frame,
Detachable desk armrest, fixed footrest,
Solid castor, solid rear wheel



CCW87

Aluminum chair frame, reclining high back, adjustable headrest, detachable armrest, elevating footrest, solid castor, PU MAG wheel, adjustable seat angle, safety belt

MAIN PARAMETERS (CM)

MODEL NO.										CTNS 20" FCL	CTNS 41" HG FCL
CCW87	52	32	40	15	50	118	75	77x40x95	20.8 23	90	223
CCW86	52	38	40	15	50	118	75	84x40x95	18.3 20.5	90	223
CCW120	43	30	30	15	48	103	100	80x41x89	16.9 18.4	72	177



CCW86

Aluminum chair frame, adjustable headrest, detachable armrest, fixed footrest, solid castor, pneumatic MAG wheel, adjustable seat angle, safety belt



CCW120

Aluminum chair frame, reclining high back, adjustable headrest, solid castor, PU MAG wheel, adjustable seat and adjustable backrest with 5-pieces seat belt

CCW88

Powder coating steel frame, reclining high back, adjustable and detachable headrest detachable pillow, adjustable seat angle height adjusted armrest, elevating footrest PU castor, PU MAG wheel, united brake



CCW89

Powder coating steel frame, reclining high back, adjustable and detachable headrest detachable pillow, adjustable seat angle, height adjusted armrest, elevating footrest PU castor, quick release pneumatic MAG wheel, united brake

MAIN PARAMETERS (CM)

MODEL NO.										CTNS 20' FCL	CTNS 40' HQ FCL
CCW88	53	46	30	20	53	118	150	82x61x84	35.5 37.7	70	148
CCW89	67	46	61	20	53	118	150	88x67x84	38.4 40.6	58	132
CCW29	60	41	60	20	50	120	90	98x35x89	23.6 26	92	230
CCW28	60	41	60	20	50	120	100	95x34x93	24.3 26.4	92	210
CCW106	56	41	40	15	50	108	100	83x33x109	20.2 22.4	98	232
CCW31	68	48	61	20	50	127	100	96x35x89	18.5 18.7	92	230



CCW29

Powder coating steel frame, reclining high back detachable armrest, elevating footrest, solid castor solid rear wheel



CCW28

Powder coating steel frame, reclining high back, movable armrest, elevating footrest, solid castor solid rear wheel



CCW106

Chromed steel frame, reclining high back, flip-up desk armrest, elevating footrest, solid castor, pneumatic MAG wheel, foldable push bar



CCW31

Aluminum chair frame, reclining high back adjustable and detachable headrest detachable armrest, elevating footrest solid castor, pneumatic rear wheel plastic hand rim

Silla de rueda
para niños



CCW64

Children type aluminum chair frame, flip-up armrest
fixed footrest, solid castor, pneumatic rear wheel
drop back handle

MAIN PARAMETERS (CM)

MODEL NO											CTNS 20' FCL	CTNS 40' HQ FCL
CCW64	51	35	57	15	48	89	75	90x33x73	13.8	16	148	315
CCW107	38	32	12.5	10	45	82	100	58x39x26	6.2	8.2	496	1200
CCW108	38	32	15	10	45	84	100	71x39x27	6	7.5	357	864
CCW109	38	32	15	12.5	45	87	100	54x20x58	7	8.5	472	864
CCW66	50	38	15	12.5	45	90	100	48x34x66	7	8.5	270	670



CCW107 ▲

Mini foldable portable aluminum chair frame
flip-up armrest, solid castor, pedal-locking
rear MAG wheel, seat belt, leg strap



CCW108 ▲

Mini foldable portable aluminum chair frame, flip-up
armrest, solid castor, pedal-locking rear MAG wheel
seat belt, leg strap



CCW109 ▲

Mini foldable portable aluminum chair frame, foldable
backrest, flip-up armrest, flip-up foot plate, solid castor,
pedal-locking PU rear wheel, seat belt, calf strap



CCW66 ▲

Transit type foldable & portable aluminum chair frame
flip-up armrest, solid castor, solid MAG wheel
locking united brake, with seat belt



CCW59

Aluminum chair frame, flip-up armrest, fixed footrest, solid castor, pneumatic rear wheel, drop back handle



CCW110

Children type, aluminum chair frame, fixed armrest, fixed footrest, solid castor, pneumatic rear wheel, drop back handle with drum brake



CCW112

Aluminum lightweight chair frame, fixed armrest, flip-up footplate, solid castor, pneumatic MAG wheel



CCW32

Aluminum chair frame, fixed armrest, solid castor, pneumatic rear wheel, fixed footrest, with united brake

MAIN PARAMETERS (CM)

MODEL NO	↔	↔	↻	↻	↻	↻	(kg)	DM	WR GR	CTNS 20' FCL	CTNS 40' HQ FCL
CCW59	57	41	61	15	50	91	100	94x33x73	15 17.2	126	284
CCW112	53	30	52	15	44	85	75	78x28x68	12.2 14.4	186	444
CCW110	62	41	58	15	48	90	100	77x30x91	11.3 13.5	108	272
CCW32	63	41	56	15	46	88	100	82x29x89	13 15	108	272
CCW111	81	41	59	20	49	87	100	103x25x91	11.7 13.7	144	372
CCW39	81	41	60	20	45	89	100	94x29x92	12.7 14.9	108	272
CCW91	81	48	58	15	45	87	100	96x27x92	11.2 13.2	108	272
CCW41	81	41	57	15	44	88	100	94x28x92	13 15	124	272



CCW111

Aluminum lightweight chair frame, fixed inclined armrest, solid castor, pneumatic rear wheel, fixed footrest



CCW39

Aluminum lightweight chair frame, fixed armrest, fixed footrest, solid castor, pneumatic rear wheel



CCW91

Aluminum lightweight chair frame, fixed armrest, fixed footrest, solid castor, pneumatic rear wheel



CCW41

Aluminum lightweight chair frame, fixed inclined armrest, fixed footrest, solid castor, pneumatic rear wheel



CCW38

Aluminum lightweight chair frame, fixed armrest, fixed footrest, solid castor, pneumatic rear wheel, plastic hand rim, drop back handle with united brake



CCW37

Aluminum lightweight chair frame, fixed armrest, fixed footrest, solid castor, pneumatic rear wheel, drop back handle, with drum brake



CCW33

Aluminum lightweight chair frame, fixed armrest, fixed footrest, solid castor, pneumatic rear wheel

MAIN PARAMETERS (CM)

MODEL NO.	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Weight	Seat Size	Seat Depth	Seat Height	Seat Width	Seat Depth
CCW38	63	41	56	15	46	88	100	94x28x75	13.2	15.4	164	300
CCW37	63	41	56	15	46	88	100	94x28x75	13.6	15.8	164	300
CCW33	63	41	56	15	46	88	100	94x28x92	13	15	108	272
CCW43	61	41	51	15	47	89	100	83x28x90	12.6	14.8	128	324
CCW65	47	37	20	12.5	45	88	100	75x25x66	6.7	8.5	228	560
CCW44	58	41	40	15	47	88	100	80x28x90	10	12	140	340
CCW45	58	41	40	15	47	89	100	80x28x90	12	15	140	340
CCW42	61	43	51	15	47	93	100	87x28x78	11	13	156	318
CCW47	56	41	55	15	50	88	100	80x28x91	13	15	140	294



CCW43

Aluminum lightweight chair frame, fixed armrest, fixed footrest, solid castor, pneumatic rear wheel with drum brake



CCW65

Mini aluminum chair frame, solid castor, solid PU MAG wheel, flip-up armrest, flip-up footplate, drop back handle, locking united brake, with seat belt



CCW44

Aluminum lightweight chair frame, fixed armrest, fixed footrest, solid castor, pneumatic MAG wheel



CCW45

Aluminum lightweight chair frame, fixed armrest, fixed footrest, solid castor, pneumatic MAG wheel, with united brake



CCW42

Super lightweight, powder coating, Pu fixed armrest, fixed footrest, solid castor, Pneumatic rear wheel, drop back handle



CCW47

Aluminum frame, flip up armrest, Detachable footrest, solid castor, Pneumatic wheel



CCW36

Aluminum lightweight chair frame, fixed armrest
fixed footrest, solid castor, pneumatic rear wheel
drop back handle, with drum brake



CCW101

Aluminum lightweight chair frame, fixed armrest
fixed footrest, solid castor, pneumatic rear MAG wheel

MAIN PARAMETERS (CM)

MODEL NO	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair
	1	2	3	4	5	6	7	8	9	10	11	12
CCW36	81	41	40	15	46	88	100	82x28x73	12	14	176	410
CCW101	81	46	30	20	50	87	100	99x26x88	9.6	11.8	108	272
CCW34	81	41	40	15	46	88	100	87x28x89	10.9	13.1	108	272
CCW80	59	46	30	15	49	90	100	79x25x91	11	13	140	324
CCW51	56	41	40	15	50	89	100	78x28x76	13.9	16.1	186	378
CCW93	58	41	40	15	50	87	100	78x28x76	13.9	16.1	186	378
CCW102	60	41	20	20	48	98	100	63x28x79	9.9	11.9	200	485
CCW53	56	41	40	15	49	90	100	84x28x90	13	15	108	272



CCW34

Aluminum lightweight chair frame, fixed armrest, fixed footrest solid castor, pneumatic rear MAG wheel, with united brake



CCW90

Aluminum lightweight chair frame, fixed armrest, flip-up footplate, solid castor, pneumatic rear MAG wheel



CCW51

Aluminum chair frame, flip-up desk armrest, detachable footrest, solid castor, pneumatic rear MAG wheel drop back handle, with drum brake



CCW93

Aluminum chair frame, flip-up inclined armrest; detachable footrest, solid castor, PU rear MAG wheel, drop back handle with united brake



CCW102

Aluminum chair frame, fixed armrest, detachable footrest solid castor, solid rear MAG wheel, drop back handle



CCW53

Aluminum chair frame, flip-up desk armrest fixed footrest, solid castor, pneumatic rear wheel



CCW92

Powder coating steel frame, flip-up armrest, detachable footrest, solid castor, solid MAG wheel



CCW46

Aluminum frame, flip-up desk armrest, detachable footrest, solid castor, pneumatic rear wheel, with united brake



CCW49

Aluminum frame, fixed armrest, detachable footrest, solid castor, PU rear wheel



CCW50

Aluminum frame, flip-up desk armrest, fixed footrest, solid castor, pneumatic rear wheel

MAIN PARAMETERS (CM)

MODEL NO	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	
	W	D	H	W	H	H	H	W	H	W	H	
CCW92	58	41	58	15	50	87	100	84x28x78	14	16.2	140	324
CCW46	56	41	56	15	50	89	100	82x28x90	13.7	15.9	108	272
CCW49	63	41	61	20	46	88	100	80x28x91	14	16	108	272
CCW50	58	41	55	15	49	85	100	94x28x80	13	15	108	272
CCW58	58	41	61	15	50	90	100	80x28x91	14.7	16.9	108	272
CCW95	64	41	61	15	51	93	100	85x38x75	17.6	20.7	135	324
CCW60	61	41	61	15	51	93	100	84x33x93	17	19	112	264
CCW61	61	41	40	15	51	91	100	78x33x93	15.7	17.7	108	272



CCW58

Aluminum frame, flip-up desk armrest, detachable footrest, solid castor, quick release pneumatic rear wheel



CCW95

Aluminum frame, drop back handle with locking drum brake, flip-up inclined armrest, height-adjustable armrest pad, elevating footrest, solid castor, quick release pneumatic rear wheel, anti wheel



CCW60

Aluminum frame, flip-up desk armrest, height-adjustable armrest pad, elevating footrest, solid castor, quick release pneumatic rear wheel, anti wheel



CCW61

Aluminum frame, flip-up desk armrest, height-adjustable armrest pad, detachable footrest, solid castor pneumatic MAG wheel, anti wheel



plaza médica
Especialistas en Equipamiento Médico

Sillas
Traumatológicas



CCW55

Aluminum chair frame, flip-up & height-adjustable desk armrest, detachable footrest, solid castor, quick release PU MAG wheel, with anti wheel, adjustable footplate



CCW62

Adjustable aluminum chair frame, flip-up inclined armrest, detachable footrest, solid castor pneumatic rear wheel



CCW56

Adjustable aluminum chair frame, flip-up & height-adjustable inclined armrest, detachable footrest, solid castor, quick release pneumatic rear wheel, with carrying wheel and anti wheel



CCW54

Aluminum chair frame, flip-up & height-adjustable inclined armrest, detachable footrest, solid castor quick release pneumatic rear wheel, with carrying wheel and anti wheel

MAIN PARAMETERS (CM)

MODEL NO	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Weight	Seat	Back	CTNS 20 FCL	CTNS 40 FCL
CCW55	61	41	61	15	52	93	100	73x33x93	17.5 19.2	92	224
CCW62	58	41	61	15	51	87	100	83x28x90	16.2 18.4	126	282
CCW56	60	41	61	15	51	87	100	84x33x93	16.2 18.4	126	282
CCW54	66	48	61	15	50	90	100	84x33x93	15.2 17.8	112	264
CCW96	67	46	60	17.5	50	93	100	76x32x94	20.4 22.6	126	282
CCW97	59	41	61	15	50	90	100	84x33x93	16.4 20.2	112	264
CCW98	68	46	60	15	50	91	100	77x30x92	17.9 20.1	108	272
CCW99	65	41	61	17.5	50	80	100	80x28x91	14 16.2	126	282



CCW96 ▲

Aluminum chair frame, flip-up & height-adjustable inclined armrest, detachable elevating footrest, PU castor quick release PU rear wheel, with anti wheel adjustable footplate



CCW97 ▲

Aluminum chair frame, flip-up & height-adjustable inclined armrest, detachable footrest, solid castor, quick release pneumatic rear wheel, with carrying wheel and anti wheel adjustable footplate, locking drum brake



CCW98 ▲

Aluminum chair frame, flip-up & height-adjustable inclined armrest, detachable footrest, PU castor, quick release PU rear wheel, locking drum brake with anti wheel, adjustable footplate



CCW99 ▲

Aluminum chair frame, flip-up & height-adjustable inclined armrest, detachable footrest, PU castor, quick release pneumatic rear wheel, with carrying wheel adjustable footplate



CCW57

Aluminum chair frame, flip up desk armrest,
Elevating footrest, solid castor,
Quick release pneumatic rear wheel



CCW206

Aluminum light weight chair frame, fixed armrest
flip-up footrest, solid castor, drop back handle
locking drum brake



CCW115

Aluminum light weight chair frame, fixed armrest, flip up footrest
solid castor, pneumatic rear MAG wheel, drop back handle
locking drum brake



CCW114

Aluminum lightweight chair frame, fixed armrest
flip-up footplate, solid castor, pneumatic rear
MAG wheel

MAIN PARAMETERS (CM)

MODEL NO	↔	↔	↻	↻	↻	↻	↻	↻	↻	↻	↻	↻
CCW57	58	41	50	15	50	90	110	80x28x91	14 16	126	260	
CCW206	55	46	30	15	50	87	100	72x30x72	11.5 13	168	384	
CCW115	80	46	61	15	48	90	100	71x28x71.5	9.5 11	192	441	
CCW114	59	46	30	15	49	90	100	72x25x76	11 13	140	324	
CCW205	70	46	58	15	49	87	100	82x28x76	11.5 13	168	401	
CCW113	64	41	61	15	51	93	100	85x38x75	17.6 20.7	135	324	
CCW204	58	46	30	15	48	94	100	72x28x72	10.6 12	192	441	
CCW116	61	46	51	15	49	92	100	82x26x76	12.3 14	168	370	



CCW205 ▲

Aluminum light weight chair frame, fixed armrest, flip-up footrest solid castor, pneumatic rear wheel, plastic hand rim, locking drum brake



CCW204 ▲

Aluminum light weight chair frame, fixed armrest fixed footrest, solid castor, pneumatic rear MAG wheel, drop back handle, locking drum brake



CCW116 ▲

Aluminum light weight chair frame, fixed armrest fixed footrest, solid castor, pneumatic rear wheel drop back handle, locking drum brake






plaza médica
Especialistas en Equipamiento Médico

Silla de ruedas



CCW69 

MAIN PARAMETERS (CM)

MODEL NO										CTNS 20' FCL	CTNS 40' HQ FCL
CCW69	67	38	61	10	46	77	100	77x61x50	9.5 - 11.8	108	230
CCW67	66	36	59	10	48	72	100	87x36x73	10.8 - 13	120	240
CCW68	60	36	59	10	48	72	100	87x36x73	11.6 - 13.6	120	240
CCW70	60	36	59	10	48	72	100	87x36x73	13 - 15	120	240
CCW118	80	36	60	10	44	77	100	87x36x79	12.4 - 14.5	120	240



CCW67



CCW118



CCW68



CCW70



CCW201



CCW202



CCW203

MAIN PARAMETERS (CM)

MODEL NO.	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	Wheelchair	CTNS X ² PCL	CTNS 4 ¹ 16 ¹ PCL
CCW201	66	36	61	10	52	81	100	78x61x57	10.1 12.1	108	230
CCW202	66	36	61	10	52	84	100	76x61x54	10 12	108	230
CCW203	60	36	61	10	48	76	100	87x36x79	11.7 13.9	120	240
CCW119	61	38	60	10	49	79	100	85x61x57	13 15	108	230



CCW73



CCW119



CCW74



CCW76

silla de rueda
deportivas



CCW72



CCW75



CCW77



CCW79




CCW80




CCW81

silla de rueda
deportivas



CCW82 



CCW83 



CCW71 



CCW78 



CCW200 

Silla de ruedas eléctrica



CPW11

Powder coating aluminum frame, detachable armrest, PU castor, pneumatic drive wheel, lithium battery, foldable chair frame

MAIN PARAMETERS (CM)

MODEL NO.	←→	←→	↻	↻	↻	↻	↻	↻	+ -	W	↻	↻	⚙️	FAV DIM	MV BW	CTNS BY PCK	CTNS BY PCK
CPW11	64	41	61	17	50	76	20	1-7	20AH	100WX2	3cm	8"	100	109x41x82	28 34.5	64	198
CPW12	60	41	61	17	50	76	20	1-7	20AH	100WX2	3cm	8"	100	106x36x78	24.5 27.8	64	198
CPW13	63	43	61	17	50	92	20	1-7	20AH	100WX2	3cm	8"	100	91x44x86	20.5 32.9	96	234
CPW14	63	43	61	17	50	92	20	1-7	20AH	100WX2	3cm	8"	100	91x44x86	20.5 32.9	96	234
CPW07	60	46	52	20	50	91	20	1-7	28AH	200WX2	4cm	8"	100	87x47x94	11.2 19.2	44	88



CPW12 ▲

Powder coating aluminum frame, detachable armrest, PU castor, pneumatic drive wheel, lithium battery, foldable chair frame



CPW13 ▲

Powder coating aluminum frame, detachable inclined armrest, detachable footrest, PU castor, pneumatic drive wheel, lithium battery, foldable chair frame



CPW14 ▲

Powder coating aluminum frame, detachable inclined armrest, detachable footrest, PU castor, pneumatic drive wheel, lithium battery, foldable chair frame



CPW07 ▲

Powder coating steel frame, detachable armrest, detachable footrest, rubber foam castor, inflat-free hollow drive wheel with hand rim, drop back handle, foldable chair frame

Silla de ruedas
eléctrica

plaza médica®
Especialistas en Equipamiento Médico

CPW02

Powder coating steel frame, drop back handle, flip-up PU desk armrest, detachable footrest, rubber foam castor and drive wheel, foldable chair frame



CPW10

Resin-coating steel frame, fixed armrest, detachable footrest, driving MAG wheel, rubber foam castor foldable chair frame

MAIN PARAMETERS (CM)

MODEL NO.	W	D	H	W	H	W	H	W	H	W	H	W	H	W	H	W	H
CPW02	63	46	31	20	50	92	20	1-6	28AH	200WX2	4cm	12"	110	83x44x85	58x63.2	50	104
CPW10	62	46	30	20	48	91	20	1-6	28AH	200WX2	4cm	12"	100	86x45x109	49 57	36	78
CPW04	63	46	31	20	53	94	20	1-6	28AH	200WX2	4cm	12"	125	83x63x74	49 57	46	96
CPW03	58	41	31	20	50	91	20	1-6	28AH	200WX2	4cm	12"	110	84x44x81	58.2 66.2	50	104
CPW15	67	43	35	22	68	100	25	1-8	38AH	250WX2	5cm	12"	120	90x74x100	75 85	28	90
CPW01	65	43	35	22	58	100	25	1-8	28AH	250WX2	5cm	12"	120	82x69x79	64.6 72	50	104



CPW04 ▲

Powder coating steel frame, flip-up desk armrest detachable footrest, rubber foam castor and drive wheel, foldable backrest



CPW03 ▲

Powder coating steel frame, drop back handle flip-up armrest, detachable footrest, rubber foam castor and drive wheel, foldable chair frame



CPW01 ▲

Powder coating steel frame, adjustable armrest pad detachable sideboard, detachable strong footrest angle-adjustable backrest, rubber foam castor and drive rear wheel, with safety lamp



CPW15 ▲

Powder coating steel frame, adjustable armrest pad detachable sideboard, detachable elevating footrest adjustable footplate, pneumatic castor and drive rear wheel, angle-adjustable backrest, with safety lamp adjustable headrest

Silla de ruedas
eléctrica

plaza médica
Especialistas en Equipamiento Médico

CPW17

Aluminum frame, drop back handle, flip-up armrest
flip-up footrest, PU castor and rear wheel, lithium
battery




CPW18

Resin-coating steel frame, adjustable and detachable
armrest PAD, detachable sideboard, rubber foam castor
and drive rear wheel, electrically lifting seat plate

MAIN PARAMETERS (CM)

MODEL NO.	W	D	H	W	H	H	W	H	+	-	W	H	H	W	D	W	D	W	D
CPW17	61	46	38	17	52	92	25	1-6	20AH	100WQ	3cm	12"	100	80x35x85	23	25	126	252	
CPW18	62	43	31	20	52	92	20	1-6	38AH	200WQ	4cm	12"	110	80x65x64	75	85	63	168	
CPW19	67	43(46)	35	23	55	113	22	1-7	38AH	200WQ	5cm	12"	110	80x70x84	82	90	42	126	



CPW19 

Powder coating steel frame, adjustable and detachable armrest PAD, detachable footrest, rubber foam castor and rear wheel, reclining backrest, with safety lights adjustable seat

Silla de ruedas eléctrica

plaza médica
Especialistas en Equipamiento Médico



CPW16

Power stand up wheelchair, powder coating frame, flip-up armrest with single controller for both driving & moving up seat. Strong rubber foam drive wheel, with comfort safety legrest belt.



CPW05


Aluminium chair frame, comfortable vehicle seat, recling high back, adjustable PU armrest, elevating footrest, rubber foam driving wheel



MAIN PARAMETERS (CM)

MODEL NO.	Seat Height	Seat Width	Seat Depth	Seat Cushion	Seat Back	Seat Backrest	Seat Backrest Height	Seat Backrest Angle	Seat Backrest Recline	Seat Backrest Recline	Seat Backrest Recline	Seat Backrest Recline	Seat Backrest Recline	Seat Backrest Recline	Seat Backrest Recline	Seat Backrest Recline	Seat Backrest Recline	Seat Backrest Recline
CPW16	83	44	31	20	57	105	25	1-8	38AH	200WX2	4cm	12°	100	118x66x80	80.4x85.2	54	144	
CPW05	65	46	31	20	33	115	20	1-6	28AH	200WX2	4cm	12°	125	84x71x81	87.1x76.3	46	96	
CPW08	60	41	31	20	50	92	20	1-6	28AH	200WX2	4cm	12°	100	72x64x85	55.63	81	216	
CPW09	60	41	31	20	46	86	20	1-6	28AH	200WX2	4cm	12°	100	72x64x85	55.63	81	216	



CPW08 

Resin-coating steel frame, adjustable armrest pad, detachable sideboard, detachable footrest, detachable backrest, rubber foam castor and drive rear wheel.



CPW09 

Resin-coating steel frame, flip-up armrest, detachable footrest, detachable backrest, rubber foam castor and drive rear wheel.

Sillas Neurológica
y Traumatológicas

plaza médica[®]
Especialistas en Equipamiento Médico



CCWR06

Foldable aluminum chair frame, reclining high back, detachable armrest, elevating footrest, solid castor, pneumatic rear wheel, soft commode seat

MAIN PARAMETERS (CM)

MODEL-NO.										CTNS 20-PCL	CTNS 10-PCL
CCWR06	64	46	61	20	52	122	100	97x41x86	17 19	90	184
CCWR03	64	44	58	20	49	87	100	94x29x89	19 21.2	120	244
CCWR01	64	44	58	20	53	129	100	107x33x78	25.8 29.3	90	184
CCWR07	64	44	58	20	55	129	100	92x34x92	25.6 28.7	90	184



CCWR03

Foldable chromed steel frame, soft commode seat, fixed armrest and fixed footrest, solid castor and solid rear wheel.



CCWR01

Foldable chromed steel frame, plastic commode seat, reclining high back, detachable armrest and elevating footrest, solid castor and rear wheel



CCWR07

Foldable chromed steel frame, plastic commode seat, reclining high back with the handle under the armrest fixed armrest and elevating footrest, solid castor and rear wheel



plaza médica
Especialistas en Equipamiento Médico

Sillas Neurológicas
y Traumatológicas



CCWR05

Aluminum chair frame, PU commode seat, flip-up desk armrest, detachable footrest, solid castor, pneumatic rear wheel



CCWR04

Foldable aluminum chair frame, soft commode seat, fixed armrest and fixed footrest, solid castor, pneumatic rear wheel



CCWR02

Foldable chromed steel frame, plastic commode seat, detachable armrest and footrest, solid castor and rear wheel




CCWR08

Foldable chromed steel frame, plastic commode seat, flip-up armrest, detachable footrest, solid castor and rear wheel

MAIN PARAMETERS (CM)


MODEL NO	↔	↔	↻	↻	↻	↻	⚖️ (kg)	📏 DIM	📏 NW / GW	📦 CTNS 20" FCL	📦 CTNS 40" HQ FCL
CCWR05	61	46	58	15	48	93	100	80x35x71	14.4 16.8	115	240
CCWR04	61	41	57	15	52	90	100	94x33x90	17 19	120	244
CCWR02	60	45	68	15	53	90	100	75x28x91	19.6 22.5	126	272
CCWR08	58	45	69	15	53	90	100	75x33x80	20.1 22.3	126	272
CCWC42	54	45	7.5	7.5	49	92	100	66x56x72	7 8.5	120	240
CCWC04	55	45	7.5	7.5	45-55	78-88	100	50x21x76	7 8.2	200	240
CCWC01	55	45	12.5	12.5	54	96	100	66x53x57	10.5 12	160	320



CCWC42 

Aluminum chair frame, fixed armrest
foldable footrest, solid rear castor
with lock, plastic commode seat



CCWC04 

Foldable chromed steel frame, fixed armrest
solid rear castor with lock, plastic commode seat



CCWC01 

Powder coating steel frame, flip down armrest
detachable footrest, solid rear castor with lock
plastic commode seat



CCWR09

Foldable chromed steel frame, soft commode seat, fixed armrest and footrest, solid castor and solid rear wheel



CCWR10

Aluminum frame, high adjustable backrest, seat frame and head rest, foldable PU armrest, angle adjustable



CCWC02

Stainless steel chair frame, detachable armrest and footrest, rear castor with lock, plastic commode seat



CCWC21

Anodize aluminum frame, fix handrest, solid castor solid rear wheel with break, plastic seat

MAIN PARAMETERS (CM)

MODEL NO												
	W	D	H	W	H	W	H	W	H	W	H	W
CCWR09	84	44	58	20	50	87	100	92x29x89	19 22	120	244	
CCWR10	55	44	12.5	12.5	5.5	112	100	78x57x71.5	15.5 20	48	104	
CCWC02	55	44	12.5	12.5	51	88	100	54x33x56	10.5 12.7	160	320	
CCWC21	60	49	7.5	10	51	90	100	64x22.5x96	8.3 10	180	378	
CCWC12	54	44	10	10	43	92-97	100	53x42x53	11 12.3	180	378	
CCWC13	55	44	12.5	12.5	51	88	100	54x53x53	12 14.2	160	320	



CCWC22

Chrome steel frame, adjustable front leg



CCWR13

Powder coating steel frame, detachable armrest and footrest, solid rear castor with lock, plastic commode seat



CCWR12

Powder coating steel frame, fix armrest, flip up footrest, solid castor, solid rear wheel with brake, plastic seat



CCWC23

MAIN PARAMETERS (CM)

MODEL NO										
CCWC22	50	43	45	81	43	33	4.1	53x14x83	1	5
CCWC23	51	41	44-60	76-92	42	33	6.4	58x23x78.5	1	7.7



CCWC03

Powder coating steel frame, fixed armrest, foldable footrest, rear castor with lock, plastic commode seat.



CCWC05

Stainless steel chair frame, PU commode seat with cover, flip down armrest, foldable footrest, detachable backrest, rear castor with lock.



CCWC07

Stainless steel chair frame, flip down armrest, detachable footrest, all four castors with lock, PU commode seat.

MAIN PARAMETERS (CM)

MODEL NO										CTNS 20' TCL	GTNS 40' HD TCL
CCWC03	55	44	12.5	12.5	47	82	100	50x26x66	10.5 12	180	378
CCWC05	54	44	50	12.5	12.5	87	100	61x55x55	12.2 14	180	368
CCWC07	54	44	12.5	12.5	51	96.5	100	54x52.5x66	11.7 13.2	145	269



CCWC14 ▲



CCWC18 ▲













CCWC41 ▲



CCWC40 ▲

MAIN PARAMETERS (CM)

MODEL NO										
	CM	CM	CM	CM	CM	kg	kg	CM	CM	kg
CCWC14	56	44	47	94	43	38	7.2	50x19x75	1	8.4
CCWC18	51	46	42-52	71-81	43	29	4.6	60x3.4x52	2	10.5
CCWC41	56	44	47	94	43	38	5.3	50x19x75	1	6.5
CCWC40	52	49	46	90	43	44	7.2	48x25.5x93	1	8.8

Comodo

plaza médica
Especialistas en Equipamiento Médico



CCWC08



CCWC43



CCWC09



CCWC10

MAIN PARAMETERS (CM)

MODEL NO.									PC CM	G.W KG
CCWC08	48	45	45-55	78-88	43	33	3.5	50x23x76	1	4.7
CCWC43	51	44	43-58	76-81	40	33	6.1	41x58x43	2	13.6
CCWC09	48	48	37-47	70-80	44	33	6	56x40x36	1	7.2
CCWC10	52	44	45-55	77-87	43	32	6.2	52x23x79	1	7.6
CCWC19	54	43	50	89	50	42	7	60x55x93	1	8
CCWC11	48	45	45-55	78-88	43	33	6	50x21x76	1	7.1
CCWC12	41	38	45	-	40	-	2.4	67x20x39	2	5.7
CCWC17	56	46	41-51	74-84	40	33	6	60x32x48	2	13.4
CCWC13	41	38	40	-	40	-	1.6	75x30x40	6	10.7
CCWC15	50	49	50-60	80-90	40	33	5	71x52x56	2	12
CCWC20	48	43	44	78	40	35	4.2	79x43x55	4	18.4
CCWC16	51	50	51-61	67-77	40	33	5.2	71x55x56	2	12.4



CCWC19



CCWC11



CCWC12



CCWC17



CCWC13



CCWC15



CCWC20



CCWC16

Comodo

plaza médica
Especialistas en Equipamiento Médico



CCWC24



CCWC25



CCWC26



CCWC27

MAIN PARAMETERS (CM)

MODEL NO.									PCB CM	G.W. (KG)
CCWC24	43	39.5	40-50	62-72	41	-	2.6	45x16x63	1	3.6
CCWC25	46	41	40-55	61-76	46	23	6.4	50x23x65	1	7.5
CCWC26	41	38	45	-	40	-	2.4	67x20x39	2	5.7
CCWC27	44	46	41-56	60-75	37	18	5.4	44x22x66	2	12
CCWC28	60	48	45	82	41	40	6.4	60.5x20x88	1	8
CCWC29	58	43	45	83	40	35	6.5	53x15.5x106.5	1	7.6
CCWC30	56	45	44	75	40	46	18.2	65x50x90	1	19.2



CCWC28



CCWC29



CCWC30



CSC11



CSC05



CSC03



CSC08



CSC12

MAIN PARAMETERS (CM)

MOBEL NO									PCB (CM)	P.W (kg)
CSC11	39	51	47	52	30	-	3.2	74x19x52	2	7.7
CSC05	41	51	39-49	76-86	37	37	4.4	80x56x40	2	11
CSC03	50	51	-	52	37	37	4.4	77x56x38.9	2	10.6
CSC08	23	45	-	18	23	-	1.5	23x19x73.5	1	2
CSC12	59	48	50	88	37	39	6.5	102x20x54	1	8
CSC07	42	50	35-45	74-84	30	40	3	68x19x39	2	6.9
CSC10	42	50	35-45	39-49	30	-	1.8	60x17x34.5	2	4.2
CSC01	50	49	49-59	82-92	37	37	4	85x55x56	2	10
CSC02	38	48	51-63	69-81	28	-	2.5	55x41x78	2	6.8
CSC13	42	50	34-45	74-84	30	40	3	68x19x39	2	7
CSC09	45	41	50-62	50-62	37	-	4.9	56x22x51	2	10.7
CSC06	51	51	34-44	73-83	37	41	3.3	53x50x25	2	7.6
CSC04	45	72	43-53	80-90	41	37	4.7	77x22x45	2	10.5



CSC07



CSC04



CSC10



CSC01



CSC02



CSC13



CSC09



CSC06

Silla de Ducha

plaza médica
Especialistas en Equipamiento Médico



CSC17



CSC15



CSC16



CSC18

MAIN PARAMETERS (CM)

MODEL NO.										
CSC14	31	66	9.5	13.5	31	-	2.8	32x7x71.5	1	3.3
CSC15	60	35	35-44	61-72	30	30	3	40x19x75	3	3.8
CSC16	42	50	42-52	75-85	30	40	3.5	40x50x43	1	4.5
CSC17	43	50	44-47.5	47.5-50.5	30	-	2.4	51x10x39	1	2.7
CSC18	42	50	42-52	45-55	30	-	2.6	40x50x43	1	3.6
CSC19	34	32.5	38-56	-	-	-	1.5	35x35x39	2	3.8
CSC20	42	63	45-55	83-93	42	38	6	71x51x91.5	2	14.5
CSC21	51	64	48-56	84-94	42	38	6.5	70x45x52	2	15
CSC22	51	41	29-14	29	23	-	3.4	50x42x46	2	7.9
CSC23	27	41	14	23	41	-	1.8	48x27x19	2	4.3



CSC14



CSC20



CSC23



CSC22



CSC12



CSC21

Caminante

plaza médica
Especialistas en Equipamiento Médico



CW26L



CW18



CW19



CW20



CW21



CW22



CW17

MAIN PARAMETERS (CM)

MODEL NO						Pcs per box	Wt. (kg)
CW26L	84	21	79-98	4.5	76x19x96	2	10.8
CW07	48	44	82-100	2.6	57x19x85	2	6.3
CW08	47	45	77-95	2.6	58x18x80	2	7
CW12	47	47	78-98	2.6	53x17x81	2	6.4
CW13	47	46	78-96	2.4	54x18x80	2	8
CW14	47	46	78-98	2.4	54x18x80	2	6
CW15	48	44	79-96	2.7	57x18x82	2	6.7



CW07



CW08



CW12



CW13



CW14



CW15



CW23



CW24



CW04



Caminante

plaza médica
Especialistas en Equipamiento Médico



CW25



CW16



CW09



CW01

MAIN PARAMETERS (CM)

MODEL NO.								
	56	44	43-60	74.5-92	2.8	57x19x79	2	8
CW26S	56	35	-	45-55	3	68x18x65	2	7.2
CW16	49	44	43-60	74.5-92	2.8	57x19x79	2	8
CW03	49	44	-	85-100	2.6	58x18x81	2	7
CW09	49	44	-	83-98	2.6	58x18x81	2	7
CW11	53	47	55-66	82-92	4	61x36x94	2	9
CW01	54	35	-	45-55	3	57x18x64	2	8
CW06	42	38	-	75-85	1.2	86x68x35	10	15



CW11



CS47



CW03



CW06



CW26S

Andador

plaza médica
Especialistas en Equipamiento Médico



MAIN PARAMETERS (CM)

MODEL NO.	WHEEL SPACING	WHEEL WIDTH	WHEEL DIAMETER	WHEEL TYPE	WHEEL MATERIAL	WHEEL COATING	WHEEL WEIGHT	WHEEL PRICE
CR01	78	46	53	78-89	6.1	62x24x82	1	7.6
CR04	72	46	60	80-92	8	60x18x80	1	9.5
CR03	71	64	61	81-93	8	60x40x78	1	9
CR05	60	46	53	78-90	5	61x19x80	1	6
CW02	46	38	-	85-95	1.2	96x68x35	10	17
CR02	59	38	-	82-97	8	63x25x70	1	9
CR06	62	46	54	80-105	8.2	62x24x82	1	9





CR05



CR06



CW31S

Length: 65cm

Weight: 18 kg

Max load: 120 kg

- adjustable -grip angle
- armrest width adjusted from 53cm-76cm
- armrest height adjusted from 71cm-110cm
- leg width adjusted from 64cm-70cm
- strong supported by gas spring
- castors with brake
- optional accessories: oxygen holder, basket, drop holder

CW31L

Length: 77.5cm

Weight: 20 kg

Max load: 200 kg

- adjustable -grip angle
- armrest width adjusted from 53cm-76cm
- armrest height adjusted from 88cm-131cm
- leg width adjusted from 70cm-85cm
- strong supported by gas spring
- castors with brake
- optional accessories: oxygen holder, basket, drop holder





CS28



CS29



CS23 CS24 CS26 CS27 CS42 CS20 CS25 CS31 CS32 CS33 CS3



CS43 CS44 CS45 CS46

MAIN PARAMETERS (CM)

MODEL NO.	HEIGHT (cm)	M.W. (kg)	FAV (mm)	ASH (cm)	G.W. (kg)
CS01	90	0.3	95x16x28	20	7
CS02	72-84	0.4	66x16x29	20	8.3
CS03	70-91	0.3	66x18x29	20	7
CS04	66-88	0.3	64x16x29	20	7
CS05	87-117.5	0.5	67x27x30	20	11
CS06	66-118.5	0.55	100x30x30	20	12
CS07	75-97	0.3	74x40x14	20	7
CS08	92-122	0.55	94x30x31	20	12
CS09	76-86	0.4	62x17x29	20	8.2
CS10	76-86	0.4	62x17x29	20	8.2
CS12	68.5-92	0.3	63x14x29	20	7
CS13	89	0.3	95x16x29	20	7
CS14	89	0.3	95x16x28	20	7
CS17	78-88	0.3	82x14x29	20	7
CS18	78-88	0.3	82x17x29	20	7
CS19	78-91	0.3	65x16x29	20	7
CS30	68.5-92	0.3	66x14x29	20	7
CS11	78-88	0.3	82x14x29	20	7
CS21	77.5-100	0.35	97x34x19	20	7
CS23	77.5-100	0.3	80x30x20	20	7
CS24	75-97	0.3	65x35x18	20	7
CS26	75-97	0.3	65x35x18	20	7
CS27	80-90	0.37	61x16x26	20	6.6
CS20	71-94	0.33	65x15x26	20	7
CS31	71-100	0.32	79x29x18	20	7.4
CS32	71-100	0.35	79x29x18	20	8
CS33	75-98	0.33	65x35x18	20	7.4
CS34	75-98	0.33	65x35x18	20	7.4
CS43	67-90	0.34	63x15x28	20	7.4
CS44	69-92	0.29	63x12x25	20	6.3
CS45	62-95	0.33	67x13x28	20	7.4
CS46	71-94	0.4	65x15x27	20	8.6







MAIN PARAMETERS (CM)

MODEL NO.	Height (cm)	Wt. (kg)	Weight (kg)	PCB (cm)	Wt. (kg)
CQC16	72-87	0.8	75x34x44	10	9
CQC01	75-97	1.1	68x28x47	10	12
CQC03	75-97	0.9	78x34x39	10	10
CQC04	75-97	0.8	78x29x43	10	9
CQC05	75-97	1.1	79x30x44	10	12
CQC06	74-96	0.9	79x30x44	10	10
CQC07	71.5-94	1.1	79x34x44	10	12
CQC08	66-89	0.85	73x31x35	10	9.5
CQC09	77-96	0.6	60x30x34	10	7
CQC10	65-87.5	0.8	73x31x35	10	9
CQC11	72.5-95	0.9	75x34x39	10	10
CQC12	72-95	1.1	79x34x43	10	12
CQC13	72-95	1.1	79x34x45	10	12
CQC14	71-94	1.1	79x34x43	10	12
CQC15	72-95	0.6	80x25x45	10	7
CQ01	134-154	0.8	141x34x70	20	20.5
CQ02	114-134	1	120x33x29	20	21.5
CS34	120	0.3	54x13x27	40	7
CQ03	127-147	1	135x34x30	20	21.5
CS39	97-117	0.8	195x34x30	20	17.5
CS40	112-143	0.75	101x33x29	20	16.5
CFA01	95-134	0.75	101x33x25	20	16.5
CS37	84	0.9	88x31x41	10	10
CS35	80	0.8	84x21x44	10	9
CS38	84	0.9	85x30x46	10	10

Scooter de movilidad y
 silla de ruedas eléctrica












CARE




CPW20 


MAIN PARAMETERS (CM)

MODEL NO.	长 x 宽 x 高										GTNS 20° FOL	GTNS 45° FOL
CPW20	112x64X113cm	15	25	6.4	6.4	33AH	200Wx2	5cm	8"	136	57	136
CPW21	115x63X115cm	25	20	6.4	6.4	38AH	200Wx2	4cm	12"	150	70	167
CPW22	91x52X80cm	20	14	4.8	4.8	17AH	200Wx2	2.5cm	8"	100	102	245
CPW23	105x64X117cm	15	15	6.4	6.4	38AH	200Wx2	4cm	8"	136	69	165



CPW21 



CPW22 




CPW23 

Scooter de movilidad y
 silla de ruedas eléctrica



CARE



CPS01 




CPS02 




CPS08 




CPS04 



CPS03 



CPS05 

MAIN PARAMETERS (CM)

MODEL NO.	型号尺寸										CTNS 20' FCL	CTNS 40' HC FCL
CPS01	105x48X90cm	20	20	6.4	6.4	12AH	180W	3cm	8"	120	87	203
CPS02	105x48X90cm	20	20	6.4	6.4	12AH	180W	3cm	8"	120	87	203
CPS03	107x48X94cm	20	20	6.4	6.4	17AH	180W	3cm	8"	120	87	203
CPS04	109x48X89cm	20	20	6.4	6.4	17AH	180W	3cm	8"	120	87	203
CPS05	125x64X117cm	25	25	8	8	38AH	400W	5cm	15"	120	43	108
CPS08	131x64X117cm	25	25	8	8	38AH	400W	5cm	15"	120	43	108

Scooter de movilidad y
silla de ruedas eléctrica



CARE



CPS06 ▲



CPS07 ▲

MAIN PARAMETERS (CM)

MODEL NO.	尺寸											
CPS06	125x64X117cm	25	25	8	8	38AH	400W	5cm	15°	120	43	108
CPS07	131x64X117cm	25	25	8	8	38AH	400W	5cm	15°	120	43	108
CPS09	131x64X117cm	25	25	8	8	38AH	400W	7cm	15°	120	43	108
CPS10	139x65X108cm	30	30	10	10	38AH	400W	7cm	15°	120	31	77
CPS11	148x69X128cm	30	30	12	12	55AH	400W	7cm	15°	205	33	82



CPS09 



CPS10 



CPS11 



plaza médica
Especialistas en Equipamiento Médico

Scooter de movilidad y
silla de ruedas eléctrica

CARE

Scooter de movilidad y
 silla de ruedas eléctrica



CARE



CPS12

MAIN PARAMETERS (CM)

MODEL NO.	型号尺寸										GTMS 20" FCL	GTMS 10" HD FCL
CPS12	156x69X133cm	36	36	14	14	75AH	1200W	8cm	15°	205	33	82
CPS13	156x69X128cm	30	30	13	13	55AH	800W	7cm	15°	205	33	82
CPS14	130x69X128cm	30	25	13	13	55AH	800W	7cm	15°	205	33	82



CPS13 



CPS14 


Silla de ruedas
rampa para scooter



CARE

Folding Ramp

- Ideal for manual/ power wheelchairs and mobility scooters
- Sturdy and strong 600lbs weight capacity
- With carrying handle for easy portability like a suitcase
- Extremely lightweight



CCR01 Folding Ramp ▲
Size: 61x74x5CM



CCR02 Folding Ramp ▲
Size: 91x74x5CM



CCR03 Folding Ramp ▲
Size: 122x74x5CM



CCR04 Folding Ramp 
Size: 152x74x5CM




CCR05 Folding Ramp 
Size: 183x74x5CM



CCR06 Folding Ramp 
Size: 181x77x10.5CM



CCR07 Folding Ramp 
Its compact design allows it to be rolled up for transport, but it can be set up in seconds. It has a durable non-slip surface and side rails that act as a bumper guide for increased safety and security.

- Three size options:
- CCR07A: 210x76.2x5CM
 - CCR07B: 150x76.2x5CM
 - CCR07C: 90x76.2x5CM

The lift chair is designed for the elderly and the leg handicapped, also for patients who have difficulty hard to stand up from the chair.

Features: ▶

- Adopts a powered lifting mechanism that lifts the entire chair up from its base, allowing the user to easily move to a standing position
- Solid hardwood frames, kiln dried and cut to precise angles to eliminate squeaks
- Double doweling at stress points and reinforced miter blocks provide extra strength and durability
- Heavy gauge tubular steel lift mechanism with 200kgs weight capacity
- Improved filler materials crafted to create comfortable seating and eliminate sagging

Three position: ▶

- Lifting, Sitting and Reclining

Motor: ▶

- German Okin motor
- Two modes: single motor or dual motors

Cover material: ▶

PU or Microfiber

Options: ▶

- Vibration massage
- Vibration massage and heating
- Kneading message



Cover color selection: ▶





CCLC01 



Zero wall
mechanism




CCLC02 

Economy
model



CCLC03 



CCLC04 



CCLC05 

Silla elevadora

plaza médica
Especialistas en Equipamiento Médico



CCLC06 



CCLC07 



CCLC08 




CCLC09 



CCLC10 



CCLC11 




CCLC12 



CCLC13 

CCLC15 




CCLC14 




CCLC16 




CCLC17 



CCLC18 





CCLC19 




CCLC20 

Color opcional



1XX: Spray painting series 2XX: Powder coating series 3XX: Aluminum anodized series 4XX: Chromed series 5: Glass bead blasting series



- | | | | | |
|-------------------------------|---------------------------------------|------------------------------------|---|------------------------------------|
| (1) 3" CASTOR | (11) ADJUSTABLE ANTI-WHEEL SET | (21) 8" PU WHEEL | (31) 20" PNEUMATIC SPOKE WHEEL SET | (41) 24" SOLID MAG WHEEL SET |
| (2) 3" SOLID CASTOR | (12) 8" WHEEL SET WITH ALUMINUM FORK | (22) 8" PNEUMATIC WHEEL | (32) 20" RUBBER PNEUMATIC MAG WHEEL SET | (42) 24" PU MAG WHEEL SET |
| (3) 4" SOLID CASTOR | (13) 8" SOLID WHEEL SET | (23) 9" SOLID WHEEL | (33) 18" PU WHEEL SET | (43) FOAM HANDLE FOR STICK |
| (4) 5" SOLID CASTOR | (14) 4" SOLID WHEEL SET | (24) 8" RUBBER FOAM WHEEL | (34) 22" SOLID SPOKE WHEEL SET | (44) FOAM HANDLE FOR GRAYON |
| (5) 6" PNEUMATIC CASTOR | (15) 4" SOLID WHEEL WITH BRAKE | (25) 12" RUBBER FOAM DRIVE WHEEL | (35) 22" 24" PLASTIC HAND RIM | (45) FOAM HANDLE FOR WALKER |
| (6) 8"x1.8" SOLID CASTOR | (16) 6" SOLID WHEEL WITH BRAKE | (26) 12" PNEUMATIC MAG WHEEL SET | (36) 24" PNEUMATIC SPOKE WHEEL SET | (46) FOAM SUPPORT FOR CRUTCH |
| (7) 8" PNEUMATIC WHEEL SET | (17) 12" SOLID WHEEL SET | (27) 12" PNEUMATIC SPOKE WHEEL SET | (37) 24" PNEUMATIC MAG WHEEL SET | (47) PU FOOTREST PAD |
| (8) 4" SOLID WHEEL SET | (18) 8" SOLID WHEEL WITH BRAKE | (28) 14" PNEUMATIC MAG WHEEL SET | (38) 24" SOLID SPOKE WHEEL SET | (48) 30x 30 (1-55-52-54) RUBBERTIP |
| (9) 6" PU WHEEL SET | (19) 8" SOLID WHEEL WITH PLASTIC FORK | (29) 16" PNEUMATIC SPOKE WHEEL SET | (39) 22" PU SPOKE WHEEL SET | (49) SHORT PU ARMREST PAD |
| (10) 8" RUBBER FOAM WHEEL SET | (20) 8" SOLID WHEEL WITH PLASTIC FORK | (30) 22" PNEUMATIC SPOKE WHEEL SET | (40) 24" PU SPOKE WHEEL SET | (50) LONG PU ARMREST PAD |

Piezas de repuesto



- | | | | | |
|---------------------------------------|---|--------------------------------------|------------------------------------|--|
| (37) INCLINED SHORT PU ARMREST PAD | (44) 45 Ø 25MM RUBBER ELEPHANT FOOTREST | (81) PLASTIC FOOTPLATE | (100) BRAKE SET WITH PLASTIC COVER | (171) COMPLETED PU PAD DESK LENGTH ARMREST |
| (38) INCLINED LONG PU ARMREST PAD | (70, 71) Ø 25MM DETACHABLE FOOTREST | (82)(83) ALUMINUM FOOTPLATE | (101) ALUMINUM LOCK SET | (172) EXTRA LENGTH STEEL ARMREST SET |
| (39) SHORT PVC ARMREST PAD | (72) Ø 22MM DETACHABLE FOOTREST | (84) RUBBER MAT FOOTPLATE | (102) STEEL LOCK SET | (173) EXTRA LENGTH ALUMINUM ARMREST SET |
| (40) LONG PVC ARMREST PAD | (73) Ø 25MM DETACHABLE FOOTREST | (85) JAPANESE TYPE PLASTIC FOOTPLATE | (103) ALUMINUM FORK | (174) SLEWING ARMREST PAD |
| (41) CANVAS LONG ARMREST PAD | (74) Ø 22MM DETACHABLE FOOTREST | (86) BEARING FOR FRONT CASTOR | (104) DRUM BRAKE WHEEL HUB | (175) TOILET SEAT WITH COVER |
| (42) PLASTIC LONG ARMREST PAD | (75, 76) ALUMINUM DETACHABLE FOOTREST | (87) BEARING FOR REAR WHEEL | (105) QUICK RELEASE HUB | (176) PLASTIC COMMODE SEAT |
| (43) DENSITY POLYETHYLENE ARMREST PAD | (77) ALUMINUM DETACHABLE FOOTREST | (88) FOLDABLE PUSH BAR | (106) Ø3 DRUM BRAKE | (177) COMMODE SEAT CUSHION |
| (44) Ø 25MM ELEPHANT FOOTREST | (78) T-SHAPE FOOTPLATE | (89) ALUMINUM LOCK SET | (107) Ø3 DRUM BRAKE | (178) COMMODE PAN |
| (45) Ø 25MM ELEPHANT FOOTREST | (79) SMALL T-SHAPE FOOTPLATE | (90) ALUMINUM LOCK SET | (108) FULL LENGTH ARMREST PAD | (179) COMMODE BUCKET |
| (46) 80Ø 22MM ELEPHANT FOOTREST | (80) JAPANESE TYPE PLASTIC FOOTPLATE | (91) STEEL HUB | (109) SLEWING ALUMINUM ARMREST PAD | |